

CENPA Antibody



PACO03355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CENPA

Uniprot

P49450

Synonyms:

CENPA; Histone H3-like centromeric protein A; Centromere autoantigen A; Centromere protein A; CENP-A

Immunogen:

Synthesized peptide derived from the N-terminal region of human CENP-A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A