

MYL12A/MYL12B/MYL9 Antibody



PACO03361

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

MYL12A/MYL12B/MYL9

Uniprot

P19105/O14950/P24844

Synonyms:

MYL12A; MLCB; MRCL3; RLC; Myosin regulatory light chain 12A; MLC-2B; Myosin RLC; Myosin regulatory light chain 2; nonsarcomeric; Myosin regulatory light chain MRCL3; MYL12B; MRCL2; MYLC2B; Myosin regulatory light chain 12B; MLC-2A; MLC-2; M

Immunogen:

Synthesized peptide derived from human MLC-2 around the non-phosphorylation site of T17/S18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A