

POLQ Antibody



PACO03362

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

POLQ

Uniprot

O75417

Synonyms:

POLQ; POLH; DNA polymerase theta; DNA polymerase eta

Immunogen:

Synthesized peptide derived from the N-terminal region of human DNA pol Theta.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A