

ZMYND11 Antibody



PACO03367

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

ZMYND11

Uniprot

Q15326

Synonyms:

ZMYND11; BS69; Zinc finger MYND domain-containing protein 11; Adenovirus 5 E1A-binding protein; Protein BS69

Immunogen:

Synthesized peptide derived from the Internal region of human BS69.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A