

DMC1 Antibody



PACO03384

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

DMC1

Uniprot

Q14565

Synonyms:

DMC1; DMC1H; LIM15; Meiotic recombination protein DMC1/LIM15 homolog

Immunogen:

Synthesized peptide derived from the Internal region of human Dmc1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A