

AKAP11 Antibody



PACO03388

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

AKAP11

Uniprot

Q9UKA4

Synonyms:

AKAP11; AKAP220; KIAA0629; A-kinase anchor protein 11; AKAP-11; A-kinase anchor protein 220 kDa; AKAP 220; hAKAP220; Protein kinase A-anchoring protein 11; PRKA11

Immunogen:

Synthesized peptide derived from the C-terminal region of human AKAP 220.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A