

# ADCY7 Antibody



PACO03392

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

ADCY7

**Uniprot**

P51828

**Synonyms:**

ADCY7; KIAA0037; Adenylate cyclase type 7; ATP pyrophosphate-lyase 7; Adenylate cyclase type VII; Adenylyl cyclase 7

**Immunogen:**

Synthesized peptide derived from the Internal region of human A Cyclase VII.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A