

SH2D1B Antibody



PACO03407

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

SH2D1B

Uniprot

O14796

Synonyms:

SH2D1B; EAT2; SH2 domain-containing protein 1B; EWS/FLI1-activated transcript 2; EAT-2

Immunogen:

Synthesized peptide derived from the Internal region of human EAT-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A