

OR56A1 Antibody



PACO03416

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

OR56A1

Uniprot

Q8NGH5

Synonyms:

OR56A1; Olfactory receptor 56A1; Olfactory receptor OR11-75

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 56A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A