

CSNK1A1L/CSNK1A1 Antibody



PACO03419

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CSNK1A1L/CSNK1A1

Uniprot

Q8N752/P48729

Synonyms:

CSNK1A1L; Casein kinase I isoform alpha-like; CKI-alpha-like; CK1; CSNK1A1; Casein kinase I isoform alpha; CKI-alpha; CK1

Immunogen:

Synthesized peptide derived from the Internal region of human Casein Kinase I alpha1/1L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A