

# MRPL52 Antibody



PACO03420

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

MRPL52

**Uniprot**

Q86TS9

**Synonyms:**

MRPL52; 39S ribosomal protein L52; mitochondrial; L52mt; MRP-L52

**Immunogen:**

Synthesized peptide derived from the Internal region of human MRP-L52.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A