

COL3A1 Antibody



PACO03434

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

COL3A1

Uniprot

P02461

Synonyms:

COL3A1; Collagen alpha-1(III) chain

Immunogen:

Synthesized peptide derived from the N-terminal region of human COL3A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A