

PACO03441

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

GAB4/GAB2

Uniprot

Q2WGN9/Q9UQC2

Synonyms:

GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like; Growth factor receptor bound protein 2-associated protein 4; GAB2; KIAA0571; GRB2-associate

Immunogen:

Synthesized peptide derived from the N-terminal region of human Gab 2/4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A