

COL5A2 Antibody



PACO03442

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

COL5A2

Uniprot

P05997

Synonyms:

COL5A2; Collagen alpha-2(V) chain

Immunogen:

Synthesized peptide derived from the N-terminal region of human COL5A2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A