

COL6A3 Antibody



PACO03444

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

COL6A3

Uniprot

P12111

Synonyms:

COL6A3; Collagen alpha-3(VI) chain

Immunogen:

Synthesized peptide derived from the Internal region of human COL6A3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A