

PACO03449

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CTTN

Uniprot

Q14247

Synonyms:

CTTN; EMS1; Src substrate cortactin; Amplaxin; Oncogene EMS1

Immunogen:

Synthesized peptide derived from human Cortactin around the non-phosphorylation site of Y466.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A