

# TPRA1 Antibody



PACO03456

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:****Gene ID:**

TPRA1

**Uniprot**

Q86W33

**Synonyms:**

TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human GPR175.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A