

OR9A2 Antibody



PACO03474

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR9A2

Uniprot

Q8NGT5

Synonyms:

OR9A2; Olfactory receptor 9A2; Olfactory receptor OR7-2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 9A2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A