

OR1N1 Antibody



PACO03498

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR1N1

Uniprot

Q8NGS0

Synonyms:

OR1N1; OR1N3; Olfactory receptor 1N1; Olfactory receptor 1-26; OR1-26; Olfactory receptor 1N3; Olfactory receptor OR9-22

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 1N1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A