

ADORA3 Antibody



PACO03505

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

ADORA3

Uniprot

P33765

Synonyms:

ADORA3; Adenosine receptor A3

Immunogen:

Synthesized peptide derived from the C-terminal region of human Adenosine A3-R.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A