

LDLRAD3 Antibody



PACO03522

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

LDLRAD3

Uniprot

Q86YD5

Synonyms:

LDLRAD3; Low-density lipoprotein receptor class A domain-containing protein 3

Immunogen:

Synthesized peptide derived from the Internal region of human LDLRAD3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A