

ABHD12B Antibody



PACO03529

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:

Gene ID:

ABHD12B

Uniprot

Q7Z5M8

Synonyms:

ABHD12B; C14orf29; Abhydrolase domain-containing protein 12B

Immunogen:

Synthesized peptide derived from the C-terminal region of human ABHD12B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A