

ABCA13 Antibody



PACO03551

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

ABCA13

Uniprot

Q86UQ4

Synonyms:

ABCA13; ATP-binding cassette sub-family A member 13

Immunogen:

Synthesized peptide derived from the Internal region of human ABCA13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A