

# ZNF600 Antibody



PACO03561

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZNF600

**Uniprot**

Q6ZNG1

**Synonyms:**

ZNF600; Zinc finger protein 600

**Immunogen:**

Synthesized peptide derived from the Internal region of human ZNF600.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A