

# DEFA3 Antibody



PACO03568

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

DEFA3

**Uniprot**

P59666

**Synonyms:**

DEFA3; DEF3; Neutrophil defensin 3; Defensin; alpha 3; HNP-3; HP-3; HP3

**Immunogen:**

Synthesized peptide derived from the Internal region of human Defensin alpha3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A