

Phospho-DOK1 (Y398) Antibody



PACO03578

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

DOK1

Uniprot

Q99704

Synonyms:

DOK1; Docking protein 1; Downstream of tyrosine kinase 1; p62(dok); pp62

Immunogen:

Synthesized peptide derived from human Dok-1 around the phosphorylation site of Y398.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A