

# Phospho-ELK1 (S383) Antibody



PACO03607

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IP

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug

**Protein Background:**

**Gene ID:**

ELK1

**Uniprot**

P19419

**Synonyms:**

ELK1; ETS domain-containing protein Elk-1

**Immunogen:**

Synthesized peptide derived from human Elk-1 around the phosphorylation site of S383.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A