

RDM1 Antibody



PACO03617

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RDM1

Uniprot

Q8NG50

Synonyms:

RDM1; RAD52B; RAD52 motif-containing protein 1; RAD52 homolog B

Immunogen:

Synthesized peptide derived from the Internal region of human RDM1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A