

EPB41 Antibody



PACO03623

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

EPB41

Uniprot

P11171

Synonyms:

EPB41; E41P; Protein 4.1; P4.1; 4.1R; Band 4.1; EPB4.1

Immunogen:

Synthesized peptide derived from human 4.1R around the non-phosphorylation site of Y660.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A