

# Phospho-AKT1/AKT2/AKT3 (Y315) Antibody



PACO03635

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

AKT1/AKT2/AKT3

**Uniprot**

P31749/P31751/Q9Y243

**Synonyms:**

AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha; AKT2; RAC-beta serine/threonine-protein kinase; Protein kinase Akt-2; Protein kinase B

**Immunogen:**

Synthesized peptide derived from human Akt around the phosphorylation site of Y315.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A