

Phospho-ERF (T526) Antibody



PACO03638

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

ERF

Uniprot

P50548

Synonyms:

ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2

Immunogen:

Synthesized peptide derived from human ERF around the phosphorylation site of T526.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A