

# ERF Antibody



PACO03639

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ERF

**Uniprot**

P50548

**Synonyms:**

ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human ERF.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A