

# Phospho-F3 (S290) Antibody



PACO03647

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

F3

**Uniprot**

P13726

**Synonyms:**

F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142

**Immunogen:**

Synthesized peptide derived from human TF around the phosphorylation site of S290.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A