

Phospho-FANCD2 (S222) Antibody



PACO03651

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

FANCD2

Uniprot

Q9BXW9

Synonyms:

FANCD2; FACD; Fanconi anemia group D2 protein; Protein FACD2

Immunogen:

Synthesized peptide derived from human FANCD2 around the phosphorylation site of S222.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A