

OR5AS1 Antibody



PACO03656

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

OR5AS1

Uniprot

Q8N127

Synonyms:

OR5AS1; Olfactory receptor 5AS1; Olfactory receptor OR11-168

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5AS1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A