

ZNF292 Antibody



PACO03680

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ZNF292

Uniprot

O60281

Synonyms:

ZNF292; KIAA0530; Zinc finger protein 292

Immunogen:

Synthesized peptide derived from the Internal region of human ZN-16.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A