

CDC42EP4 Antibody



PACO03721

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CDC42EP4

Uniprot

Q9H3Q1

Synonyms:

CDC42EP4; BORG4; CEP4; Cdc42 effector protein 4; Binder of Rho GTPases 4

Immunogen:

Synthesized peptide derived from the Internal region of human Cdc42EP4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A