

SMUG1 Antibody



PACO03722

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

SMUG1

Uniprot

Q53HV7

Synonyms:

SMUG1; Single-strand selective monofunctional uracil DNA glycosylase

Immunogen:

Synthesized peptide derived from the Internal region of human SMUG1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A