

Phospho-DAPK2 (S318) Antibody



PACO03728

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

DAPK2

Uniprot

Q9UIK4

Synonyms:

DAPK2; Death-associated protein kinase 2; DAP kinase 2; DAP-kinase-related protein 1; DRP-1

Immunogen:

Synthesized peptide derived from human DAPK2 around the phosphorylation site of S318.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A