

Phospho-ALOX5 (S523) Antibody



PACO03746

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ALOX5

Uniprot

P09917

Synonyms:

ALOX5; LOG5; Arachidonate 5-lipoxygenase; 5-LO; 5-lipoxygenase

Immunogen:

Synthesized peptide derived from human 5-LO around the phosphorylation site of S523.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A