

RAD54B Antibody



PACO03788

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

RAD54B

Uniprot

Q9Y620

Synonyms:

RAD54B; DNA repair and recombination protein RAD54B; RAD54 homolog B

Immunogen:

Synthesized peptide derived from the Internal region of human Rad54B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A