

ABHD14A Antibody



PACO03792

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ABHD14A

Uniprot

Q9BUJ0

Synonyms:

ABHD14A; Abhydrolase domain-containing protein 14A

Immunogen:

Synthesized peptide derived from the C-terminal region of human ABHD14A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A