

# TAS2R45 Antibody



PACO03800

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TAS2R45

**Uniprot**

P59539

**Synonyms:**

TAS2R45; GPR59; Taste receptor type 2 member 45; T2R45; G-protein coupled receptor 59

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human T2R45.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A