

GPR124 Antibody



PACO03806

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ADGRA2

Uniprot

Q96PE1

Synonyms:

GPR124; KIAA1531; TEM5; G-protein coupled receptor 124; Tumor endothelial marker 5

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR124.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A