

PACO03810

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

FDCSP

Uniprot

Q8NFU4

Synonyms:

FDCSP; C4orf7; Follicular dendritic cell secreted peptide; FDC secreted protein; FDC-SP

Immunogen:

Synthesized peptide derived from the C-terminal region of human FDC-SP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A