

# Phospho-GATA1 (S142) Antibody



PACO03825

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IP

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug

**Protein Background:****Gene ID:**

GATA1

**Uniprot**

P15976

**Synonyms:**

GATA1; ERYF1; GF1; Erythroid transcription factor; Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein

**Immunogen:**

Synthesized peptide derived from human GATA-1 around the phosphorylation site of S142.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A