

Phospho-GATA3 (S308) Antibody



PACO03828

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

GATA3

Uniprot

P23771

Synonyms:

GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3

Immunogen:

Synthesized peptide derived from human GATA-3 around the phosphorylation site of S308.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A