

# NKX2-8 Antibody



PACO03829

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

NKX2-8

**Uniprot**

Q8IUT7

**Synonyms:**

Homeobox protein Nkx-2.8

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Nkx-2.8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A