

GATA6 Antibody



PACO03832

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

GATA6

Uniprot

Q92908

Synonyms:

GATA6; Transcription factor GATA-6; GATA-binding factor 6

Immunogen:

Synthesized peptide derived from the C-terminal region of human GATA-6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A