
Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

GCKR

Uniprot

Q14397

Synonyms:

GCKR; Glucokinase regulatory protein; Glucokinase regulator

Immunogen:

Synthesized peptide derived from the Internal region of human GCKR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A